

# Paediatric Normal Values – Anaesthesia

All codes listed are Intersurgical product codes

Interact with us  
  
[www.intersurgical.co.uk](http://www.intersurgical.co.uk)



Age	Weight (kg)	Height or length (cm)	Anaesthetic facemask size	Oropharyngeal airway size (ISO) *	i-gel size	LM size	LM cuff maximum inflation value (ml)	Tracheal tube uncuffed internal diameter (mm) **	Tracheal tube cuffed internal diameter (mm)	Awake heart rate (per minute)	Systolic blood pressure (mm Hg)	Respiratory rate (per minute)	Estimated tidal volume (ml)	Approximate blood volume (ml)	Acceptable haematocrit	Intravenous maintenance fluid (ml/hr) ***	Breathing system HMEF	Anaesthetic circuit
< 1 month	3.5	50	0 7290000	00 (5.0) 1110050	1 8201000	1 8001000	< 4	3.0		85 - 205	60	30 - 40	21	315	≥ 0.30	14	Clear-Therm Micro (1441000)	Jackson-Rees T-piece with 0.5 litre reservoir bag (2121000, 2122000 - with APL valve) in anaesthetic room and 15 mm paediatric circle system (2142000) in theatre
1 month	4	54						3.0 - 3.5	3.0		80		24			360		
3 months	5	60	1 7291000	0 (5.5) 1110055	1.5 8215000	1.5 8015000	< 7	3.5	3.0	100 - 180	80	30 - 40	30	400	≥ 0.25	20	Clear-Therm Micro (1441000)	
6 months	7	67						4.0	3.5		80		42			560		
1 year	10	78	2 7292000 7292001	1 (6.5) 1111065	1.5	2 8002000	< 10	4.0	3.5	100 - 180	92	24 - 30	60	800	≥ 0.2	40	Clear-Therm Mini (1831000)	
2 years	12	87			4.5			4.0	94		72		900			44		
3 years	14	95			4.5 - 5.0			4.0 - 4.5	96		84		1050			48		
4 years	16	103			5.0			4.5	98		96		1200			52		
5 years	18	109			5.0 - 5.5			4.5 - 5.0	100		108		1350			56		
6 years	20	116	3 7293000 7293001	1.5 (7.0) 1111570	2	2.5 8025000	< 14	5.5	5.0	60 - 140	102	20 - 24	120	1500	≥ 0.2	60	Clear-Therm Mini (1831000)	
7 years	22	122			5.5 - 6.0			5.0 - 5.5	104		132		1650			62		
8 years	26	128			6.0 - 6.5				106		156		1820			66		
10 years	30	139	3 8003000	2.5	2.5	3 8003000	< 20	7.0	7.0	60 - 100	110	20 - 24	180	2100	≥ 0.2	70	Clear-Therm 3 (1541000)	
12 years	38	149			7.0 - 7.5				114		228		2660			78		
Adolescent	50	161	2 (8.0) 1112080	3 8203000	4 8004000	< 30	7 - 8				118	12 - 20	300	3500		90		
Reference	1	2	3	3	3	3	3	1	1	1	1	1	1	4	4	6	3	3

Ambient temperature to be a minimum of 21 degrees Celsius. For smaller children and neonates undergoing surgery or resuscitation, additional warming with a Bair Hugger® or similar device is recommended (5).  
 \* The correct size of an oropharyngeal airway is one that, when laid against the side of the face, has a length equal to the distance between the level of the patient's incisors (or where they will be) to the angle of the jaw (1).  
 \*\* A correctly sized uncuffed tracheal tube should have a small audible leak around the tube when 20 cm of water pressure is applied from the breathing system (5).

### \*\*\* Intravenous maintenance fluid recommendations for previously well children aged from one month to 16 years old

The majority of children may be safely administered sodium chloride 0.45% with glucose (2.5 or 5%). Do not use sodium chloride 0.18% with glucose 4%.  
 Some children at high risk of hyponatraemia should only receive isotonic fluids (see list opposite).  
 Some acutely ill children with increased anti-diuretic hormone (ADH) secretion (e.g. post-operative patients or those with intracranial infections or head injuries) may benefit from their maintenance fluid being restricted to two-thirds normal recommended volume.  
 To avoid dangerous hypo or hypernatraemia, monitor the child's weight and calculate fluid balance. Use a volumetric pump. Check plasma electrolyte and glucose concentration before and regularly throughout intravenous therapy.  
 Consider adding potassium 40 mmol/l to maintenance fluids once plasma potassium levels are known.  
 Children requiring both maintenance fluids and replacement of ongoing losses should receive a single isotonic fluid.

### Children who should only receive isotonic fluids include those who:

are peri- or post-operative  
 have low plasma sodium  
 have CNS infection or a head injury  
 have sepsis  
 have a self-wasting syndromes  
 require the replacement of ongoing losses  
 have intravascular volume depletion or hypotension  
 have bronchiolitis  
 have excessive gastrointestinal losses  
 have a chronic condition such as diabetes, cystic fibrosis or a pituitary deficit  
**Examples of isotonic fluids are: sodium chloride 0.9%, sodium chloride 0.9% with 5% glucose or Hartmann's solution.**  
**For further information regarding the treatment of shock and the replacement of pre-existing fluid deficit, consult the NPSA website, EPALS manual and other appropriate resources.**

### References

- European Paediatric Life Support, 4th Edition, 2016; p3-4, 37, 149. Reproduced with the kind permission of the Resuscitation Council (UK).
- Child Growth Foundation (Charity Reg No 274325). Boys and Girls Growth Charts (Birth-18 years). London: 2 Mayfield Avenue, London. W4 1PW, 1996.
- Intersurgical product information <http://www.intersurgical.co.uk/>. Accessed June 2017.
- Cunliffe M. Fluid and electrolyte management in children. BJA CEPD reviews 2003; 3(1): p1-4.
- Basic techniques for anaesthesia. In Sumner E and Hatch DJ, eds. Paediatric Anaesthesia. London: Arnold, a member of the Hodder Headline Group 2000; p182, 194.
- Reducing the risk of hyponatraemia when administering intravenous infusions to children, March 2007; p11. <http://www.nrls.npsa.nhs.uk/resources/?EntryId45=59809>. Accessed June 2017.

### Disclaimer

Whilst every care has been taken to ensure that doses and recommendations are correct, the responsibility for final checking must rest with the practitioner. The authors cannot accept any responsibility for errors in this publication. Equipment sizes are based upon Intersurgical recommendations.

### Copyright

© Intersurgical 2018, all rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior written permission of the copyright holder.

### Accreditation

Intersurgical would like to thank Richard Hixson, Consultant Anaesthetist at Darlington Memorial Hospital for the concept and his help in the production of this chart.