



PEEP valves and accessories for the StarMed hoods



Critical Care ▪ Patient Interfaces ▪ **STARMed**

PEEP valves

Adjustable PEEP valves for CPAP therapy



Code	Connection	Box Qty.
SPV 20	22mm male	10
SPV22F	22mm female	10
SPV30F	30mm female	10

Over pressure valve

Over pressure valve for



Code	Description	Box Qty.
CP922	Over pressure valve	10

Accessories



STR/028

STR/014

STR/019

CP903

CP908

CP917

SPV108

CP924

CP900

CP920

CP906

CP912

Code	Description	Box Qty.
STR/028	Nebuliser accessory, complete with mouthpiece for delivering nebulised drugs during CPAP or NIV therapy	5
STR/014	22F/22F straight adaptor for connecting the hood to some breathing systems configurations	10
STR/019	22F cap to occlude the expiratory side of the hood in some breathing systems configurations	10
CP903	Cuff inflation bulb for inflating the built in neck cushion present in some hood ranges	5
CP908	Earplugs for improving patient comfort during therapy	25
CP917	22F/22F adaptor with built in manometer to convert any NIV hood into a CPAP hood	5
SPV108	22F/22F adaptor with built in manometer and adjustable PEEP valve to convert any NIV hood for immediate use in CPAP therapy	10
CP924	VentuPlus Venturi flow driver accessory	10
CP900	Spare underarm straps. Spare standard fixations for adult CaStar and Ventukit hoods	5
CP920	Spare underarm straps for CaStar Up. Spare fixations for adult CaStar Up openable hood series	5
CP906	Elastic belts with slots for bed fixation. An alternative hood fixation system	5
CP912	Infant harness. Spare fixation for the Infant Low hood. Available also for Infant High (ref. CP913)	5



www.intersurgical.co.uk/info/starmed

Make an enquiry



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: 0118 965 6300 F: 0118 965 6356 info@intersurgical.co.uk www.intersurgical.co.uk



Manufacturer:
Intersurgical S.p.A.
Mirandola (MO) Italy

Please think before you print
Save energy and paper. If you must print this information sheet please print it double sided.

IS1.47_UK • Issue 4 06.20