


Intersurgical InTube®

A range of ET tubes for adult and paediatric use



Introducing the new Intersurgical InTube® range of ET tubes

The Intersurgical InTube range of endotracheal tubes offer three varieties of ET tube: cuffed, uncuffed and reinforced in a comprehensive range of sizes from paediatric to adult. These are designed for use in anaesthesia, A&E, intensive care and emergency medicine.

High volume, low pressure (HVLP) cuffs

Intersurgical InTube cuffs are high volume, low pressure, providing even cuff inflation, reduced trauma and stabilisation of the tube in the centre of the trachea.

Features and Benefits



Closely spaced wire coil (on wire-reinforced only)


The reinforced range offers a closely spaced wire coil to increase flexibility and to reduce the risk of biting down

Highest quality valve

To ensure continual cuff integrity

Soft tip

To reduce unnecessary trauma



High quality 15mm connector

Reliable connection to all standard equipment

Precise calibration

Reliably indicating depth of insertion

Radio-opaque line

To assist in radiographic imaging

Murphy's eye

To reduce the risk of occlusion and to help maintain gaseous flow

Ordering information – all Intersurgical InTube ET tubes come in a box quantity of 10

Cuffed size	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0	9.5	10.0
Code	8040045	8040050	8040055	8040060	8040065	8040070	8040075	8040080	8040085	8040090	8040095	8040100

Uncuffed size	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0
Code	8050025	8050030	8050035	8050040	8050045	8050050	8050055	8050060	8050065	8050070

Reinforced size	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0	9.5	10.0
Code	8060045	8060050	8060055	8060060	8060065	8060070	8060075	8060080	8060085	8060090	8060095	8060100

IS6.4 • Issue 2 11.09